


Cuspidiol

Chemical Properties

CAS No. : 51593-96-5
Formula: C₁₄H₂₀O₃
Molecular Weight: 236.31
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Cuspidiol shows cytotoxicity against Jurkat T cell clone E6.1 at 7.3 ug/mL.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2317 mL	21.1586 mL	42.3173 mL
5 mM	0.8463 mL	4.2317 mL	8.4635 mL
10 mM	0.4232 mL	2.1159 mL	4.2317 mL
50 mM	0.0846 mL	0.4232 mL	0.8463 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481