


## Dodoviscin J

### Chemical Properties

|                   |                                                          |
|-------------------|----------------------------------------------------------|
| CAS No. :         | 1372527-42-8                                             |
| Formula:          | C <sub>22</sub> H <sub>22</sub> O <sub>7</sub>           |
| Molecular Weight: | 398.41                                                   |
| Appearance:       | no data available                                        |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year |



### Biological Description

|               |                                                                                                                      |
|---------------|----------------------------------------------------------------------------------------------------------------------|
| Description   | Dodoviscin J can promote adipocyte differentiation as characterized by increased triglyceride levels in 3T3L1 cells. |
| Targets(IC50) | transporter                                                                                                          |

### Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 2.510 mL  | 12.5499 mL | 25.0998 mL |
| 5 mM  | 0.502 mL  | 2.510 mL   | 5.020 mL   |
| 10 mM | 0.251 mL  | 1.255 mL   | 2.510 mL   |
| 50 mM | 0.0502 mL | 0.251 mL   | 0.502 mL   |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481