


F1839-I

Chemical Properties

CAS No. : 159096-49-8
Formula: C₂₃H₃₂O₄
Molecular Weight: 372.5
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	F1839-I is a natural product from <i>Stachybotrys</i> sp.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6846 mL	13.4228 mL	26.8456 mL
5 mM	0.5369 mL	2.6846 mL	5.3691 mL
10 mM	0.2685 mL	1.3423 mL	2.6846 mL
50 mM	0.0537 mL	0.2685 mL	0.5369 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481