

Gynuramide II

Chemical Properties


CAS No. : 295803-03-1

Formula: C42H83NO5

Molecular Weight: 682.11

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Gynuramide II is a natural product of Exochorda, Rosaceae. The catalog number is TN4176 and the CAS number is 295803-03-1. Gynuramide II can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.466 mL	7.3302 mL	14.6604 mL
5 mM	0.2932 mL	1.466 mL	2.9321 mL
10 mM	0.1466 mL	0.733 mL	1.466 mL
50 mM	0.0293 mL	0.1466 mL	0.2932 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481