Data Sheet (Cat.No.TN4189)



Hedychenone

Chemical Properties


CAS No.: 56324-54-0

Formula: C20H26O2

Molecular Weight: 298.42

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description Hedychenone has anti-inflammatory activity, it also shows potent in vitro cytotoxic

activity against cancerous cells.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.351 mL	16.7549 mL	33.5098 mL
5 mM	0.6702 mL	3.351 mL	6.702 mL
10 mM	0.3351 mL	1.6755 m Ĺ	3.351 mL
50 mM	0.067 mL	0.3351 mL	0.6702 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com