


Kadsurenin D

Chemical Properties

CAS No. : 140669-89-2
Formula: C₂₁H₂₄O₅
Molecular Weight: 356.41
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Kadsurenin D is a natural product for research related to life sciences. The catalog number is TN4365 and the CAS number is 140669-89-2.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8058 mL	14.0288 mL	28.0576 mL
5 mM	0.5612 mL	2.8058 mL	5.6115 mL
10 mM	0.2806 mL	1.4029 mL	2.8058 mL
50 mM	0.0561 mL	0.2806 mL	0.5612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481