


Lucidenic acid L

Chemical Properties

CAS No. : 110267-45-3
Formula: C₂₇H₃₈O₇
Molecular Weight: 474.59
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Lucidenic acid L is a natural product for research related to life sciences. The catalog number is TN4458 and the CAS number is 110267-45-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1071 mL	10.5354 mL	21.0708 mL
5 mM	0.4214 mL	2.1071 mL	4.2142 mL
10 mM	0.2107 mL	1.0535 mL	2.1071 mL
50 mM	0.0421 mL	0.2107 mL	0.4214 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481