


Paxiphylline D

Chemical Properties

CAS No. :	1092555-02-6
Formula:	C ₂₃ H ₂₉ NO ₄
Molecular Weight:	383.48
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Paxiphylline D is a natural product of Daphniphyllum, Daphniphyllaceae. The catalog number is TN4744 and the CAS number is 1092555-02-6. Paxiphylline D can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6077 mL	13.0385 mL	26.077 mL
5 mM	0.5215 mL	2.6077 mL	5.2154 mL
10 mM	0.2608 mL	1.3038 mL	2.6077 mL
50 mM	0.0522 mL	0.2608 mL	0.5215 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481