


Peucedanocoumarin II

Chemical Properties

CAS No. : 130464-56-1
Formula: C₂₁H₂₂O₇
Molecular Weight: 386.4
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Peucedanocoumarin II shows only partial relaxation.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.588 mL	12.940 mL	25.8799 mL
5 mM	0.5176 mL	2.588 mL	5.176 mL
10 mM	0.2588 mL	1.294 mL	2.588 mL
50 mM	0.0518 mL	0.2588 mL	0.5176 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481