


Pterokaurene L3

Chemical Properties

CAS No. : 77658-38-9
Formula: C₂₀H₃₀O₃
Molecular Weight: 318.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pterokaurene L3 is a natural product for research related to life sciences. The catalog number is TN4854 and the CAS number is 77658-38-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1402 mL	15.7011 mL	31.4021 mL
5 mM	0.628 mL	3.1402 mL	6.2804 mL
10 mM	0.314 mL	1.5701 mL	3.1402 mL
50 mM	0.0628 mL	0.314 mL	0.628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481