


Qianhu coumarin G

Chemical Properties

CAS No. : 68692-61-5
Formula: C₁₄H₁₄O₅
Molecular Weight: 262.26
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Qianhu coumarin G is a natural product from <i>Peucedanum praeruptorum</i> .
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.813 mL	19.065 mL	38.1301 mL
5 mM	0.7626 mL	3.813 mL	7.626 mL
10 mM	0.3813 mL	1.9065 mL	3.813 mL
50 mM	0.0763 mL	0.3813 mL	0.7626 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481