


Securinol A

Chemical Properties

CAS No. : 5008-48-0
Formula: C₁₃H₁₇NO₃
Molecular Weight: 235.28
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Securinol A is a natural product from <i>Securinega suffruticosa</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2503 mL	21.2513 mL	42.5026 mL
5 mM	0.8501 mL	4.2503 mL	8.5005 mL
10 mM	0.425 mL	2.1251 mL	4.2503 mL
50 mM	0.085 mL	0.425 mL	0.8501 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481