


Stigmastane-3,5,6-triol

Chemical Properties

CAS No. :	20835-91-0
Formula:	C ₂₉ H ₅₂ O ₃
Molecular Weight:	448.72
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Stigmastane-3,5,6-triol and derivatives could be used for preparing antiviral medicine.
Targets(IC50)	Antifection

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2286 mL	11.1428 mL	22.2856 mL
5 mM	0.4457 mL	2.2286 mL	4.4571 mL
10 mM	0.2229 mL	1.1143 mL	2.2286 mL
50 mM	0.0446 mL	0.2229 mL	0.4457 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481