


Withaperuvn C

Chemical Properties

CAS No. : 81644-34-0
Formula: C₂₈H₃₈O₇
Molecular Weight: 486.6
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Withaperuvn C may have anti-inflammatory activity with lipopolysaccharide (LPS)-activated murine macrophage RAW 264.7 cells.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0551 mL	10.2754 mL	20.5508 mL
5 mM	0.411 mL	2.0551 mL	4.1102 mL
10 mM	0.2055 mL	1.0275 mL	2.0551 mL
50 mM	0.0411 mL	0.2055 mL	0.411 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481