


Epimedonin I

Chemical Properties

CAS No. : 222285-83-6
Formula: C₂₁H₂₂O₈
Molecular Weight: 402.39
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Epimedonin I is a natural product for research related to life sciences. The catalog number is TN5662 and the CAS number is 222285-83-6.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4852 mL	12.4258 mL	24.8515 mL
5 mM	0.497 mL	2.4852 mL	4.9703 mL
10 mM	0.2485 mL	1.2426 mL	2.4852 mL
50 mM	0.0497 mL	0.2485 mL	0.497 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481