


Ganomycin I

Chemical Properties

CAS No. : 1191255-15-8
Formula: C₂₁H₂₆O₄
Molecular Weight: 342.435
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ganomycin I inhibits the activity of HIV-1 protease.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9202 mL	14.6011 mL	29.2022 mL
5 mM	0.584 mL	2.9202 mL	5.8404 mL
10 mM	0.292 mL	1.4601 mL	2.9202 mL
50 mM	0.0584 mL	0.292 mL	0.584 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481