


Bombiprenone

Chemical Properties

CAS No. : 21978-49-4
Formula: C43H70O
Molecular Weight: 603.032
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Bombiprenone is a natural product of Momordica, Cucurbitaceae. The catalog number is TN5849 and the CAS number is 21978-49-4. Bombiprenone can be used as a reference standard.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6583 mL	8.2915 mL	16.5829 mL
5 mM	0.3317 mL	1.6583 mL	3.3166 mL
10 mM	0.1658 mL	0.8291 mL	1.6583 mL
50 mM	0.0332 mL	0.1658 mL	0.3317 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481