


Lipodeoxyaconitine

Chemical Properties

CAS No. : 244190-83-8
Formula: C₅₀H₇₅NO₁₀
Molecular Weight: 850.147
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Lipodeoxyaconitine is a natural product for research related to life sciences. The catalog number is TN6114 and the CAS number is 244190-83-8.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1763 mL	5.8813 mL	11.7626 mL
5 mM	0.2353 mL	1.1763 mL	2.3525 mL
10 mM	0.1176 mL	0.5881 mL	1.1763 mL
50 mM	0.0235 mL	0.1176 mL	0.2353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481