


Hapepunine

Chemical Properties

CAS No. : 68422-01-5
Formula: C₂₈H₄₇NO₂
Molecular Weight: 429.68
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Hapepunine is a natural product.
-------------	----------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3273 mL	11.6366 mL	23.2731 mL
5 mM	0.4655 mL	2.3273 mL	4.6546 mL
10 mM	0.2327 mL	1.1637 mL	2.3273 mL
50 mM	0.0465 mL	0.2327 mL	0.4655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481