


Neokurarinol

Chemical Properties

CAS No. : 52483-00-8
Formula: C₂₇H₃₄O₇
Molecular Weight: 470.55
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Neokurarinol is a natural product.
-------------	------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1252 mL	10.6259 mL	21.2517 mL
5 mM	0.425 mL	2.1252 mL	4.2503 mL
10 mM	0.2125 mL	1.0626 mL	2.1252 mL
50 mM	0.0425 mL	0.2125 mL	0.425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481