


(-)-Isobicyclogermacrenal

Chemical Properties

CAS No. : 73256-82-3
Formula: C₁₅H₂₂O
Molecular Weight: 218.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(-)-Isobicyclogermacrenal is a natural product.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.580 mL	22.9001 mL	45.8001 mL
5 mM	0.916 mL	4.580 mL	9.160 mL
10 mM	0.458 mL	2.290 mL	4.580 mL
50 mM	0.0916 mL	0.458 mL	0.916 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481