


(2S)-1-O-p-coumaroyl-2-O-acetyl-3-O-β-D-glucopyranosylglycerol

Chemical Properties

CAS No. :
Formula: C₂₂H₃₀O₁₁
Molecular Weight: 470.47
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(2S)-1-O-p-coumaroyl-2-O-acetyl-3-O-β-D-glucopyranosylglycerol is a useful research compound for many research applications, including the fields of chemistry and life sciences, among others.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 50 mg/mL (106.28 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1255 mL	10.6277 mL	21.2553 mL
5 mM	0.4251 mL	2.1255 mL	4.2511 mL
10 mM	0.2126 mL	1.0628 mL	2.1255 mL
50 mM	0.0425 mL	0.2126 mL	0.4251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481