


2-Iodoadenosine

Chemical Properties

CAS No. : 35109-88-7
Formula: C₁₀H₁₂IN₅O₄
Molecular Weight: 393.14
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Nucleosides Derivatives -Halo-nucleosides; Scaffolds and Templates
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5436 mL	12.7181 mL	25.4362 mL
5 mM	0.5087 mL	2.5436 mL	5.0872 mL
10 mM	0.2544 mL	1.2718 mL	2.5436 mL
50 mM	0.0509 mL	0.2544 mL	0.5087 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481