


N4-Desmethyl wyosine

Chemical Properties

CAS No. : 59327-60-5
Formula: C₁₃H₁₅N₅O₅
Molecular Weight: 321.29
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	N4-Desmethyl wyosine is a Nucleoside Derivative - Bi-/Tri-cyclic nucleoside; naturally modified ribo-nucleoside.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1125 mL	15.5623 mL	31.1245 mL
5 mM	0.6225 mL	3.1125 mL	6.2249 mL
10 mM	0.3112 mL	1.5562 mL	3.1125 mL
50 mM	0.0622 mL	0.3112 mL	0.6225 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481