


4'-C-Methyl-5-methoxyuridine

Chemical Properties

CAS No. : 2305415-71-6
Formula: C₁₁H₁₆N₂O₇
Molecular Weight: 288.25
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside Derivatives - 4'-Modified nucleosides; 5-Modified pyrimidine nucleosides
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4692 mL	17.3461 mL	34.6921 mL
5 mM	0.6938 mL	3.4692 mL	6.9384 mL
10 mM	0.3469 mL	1.7346 mL	3.4692 mL
50 mM	0.0694 mL	0.3469 mL	0.6938 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481