


5'-O-(4,4'-Dimethoxytrityl)-3'-O-t-butyl dimethylsilyl adenosine

Chemical Properties

CAS No. :	81794-12-9
Formula:	C ₃₇ H ₄₅ N ₅ O ₆ Si
Molecular Weight:	683.87
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside; Used for special nucleoside or RNA modification
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4623 mL	7.3113 mL	14.6227 mL
5 mM	0.2925 mL	1.4623 mL	2.9245 mL
10 mM	0.1462 mL	0.7311 mL	1.4623 mL
50 mM	0.0292 mL	0.1462 mL	0.2925 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481