


6-Chloro-9-(tetrahydro-2-pyranyl)purine

Chemical Properties

CAS No. : 7306-68-5
Formula: C₁₀H₁₁ClN₄O
Molecular Weight: 238.67
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Intermediate and building block-electrophile, heterocyclic compound-purine
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1899 mL	20.9494 mL	41.8989 mL
5 mM	0.838 mL	4.1899 mL	8.3798 mL
10 mM	0.419 mL	2.0949 mL	4.1899 mL
50 mM	0.0838 mL	0.419 mL	0.838 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481