


5'-O-(4,4'-Dimethoxytrityl)-N6-Phenoxyacetyl adenosine

Chemical Properties

CAS No. :	121076-16-2
Formula:	C39H37N5O8
Molecular Weight:	703.74
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	5'-O-(4,4'-Dimethoxytrityl)-N6-Phenoxyacetyl adenosine is a 2/6/8-modified purine nucleoside.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.421 mL	7.1049 mL	14.2098 mL
5 mM	0.2842 mL	1.421 mL	2.842 mL
10 mM	0.1421 mL	0.7105 mL	1.421 mL
50 mM	0.0284 mL	0.1421 mL	0.2842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481