


2-Methoxy-2'-deoxyadenosine

Chemical Properties

CAS No. : 24757-70-8
Formula: C₁₁H₁₅N₅O₄
Molecular Weight: 281.27
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	2-Methoxy-2'-deoxyadenosine is a Nucleoside Derivative - 2-Modified purine nucleoside.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5553 mL	17.7765 mL	35.553 mL
5 mM	0.7111 mL	3.5553 mL	7.1106 mL
10 mM	0.3555 mL	1.7777 mL	3.5553 mL
50 mM	0.0711 mL	0.3555 mL	0.7111 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481