


## 2',3',5'-Tri-O-acetyl-5-cyanouridine

### Chemical Properties

CAS No. : 1630743-99-5  
Formula: C<sub>16</sub>H<sub>17</sub>N<sub>3</sub>O<sub>9</sub>  
Molecular Weight: 395.32  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



### Biological Description

|             |                                                                                                     |
|-------------|-----------------------------------------------------------------------------------------------------|
| Description | 2',3',5'-Tri-O-acetyl-5-cyanouridine is a Nucleoside Derivative - 5-Modified pyrimidine nucleoside. |
|-------------|-----------------------------------------------------------------------------------------------------|

### Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 2.5296 mL | 12.648 mL | 25.296 mL |
| 5 mM  | 0.5059 mL | 2.5296 mL | 5.0592 mL |
| 10 mM | 0.253 mL  | 1.2648 mL | 2.5296 mL |
| 50 mM | 0.0506 mL | 0.253 mL  | 0.5059 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481