


2',3'-Dideoxyuridine

Chemical Properties

CAS No. : 5983-09-5
Formula: C₉H₁₂N₂O₄
Molecular Weight: 212.2
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside Derivatives - 2',3'-Dideoxy nucleoside
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7125 mL	23.5627 mL	47.1254 mL
5 mM	0.9425 mL	4.7125 mL	9.4251 mL
10 mM	0.4713 mL	2.3563 mL	4.7125 mL
50 mM	0.0943 mL	0.4713 mL	0.9425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481