


2'-Deoxy-N4-methylcytidine

Chemical Properties

CAS No. :	22882-02-6
Formula:	C10H15N3O4
Molecular Weight:	241.24
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	2'-Deoxy-N4-methylcytidine (N(3)-Methyl-2'-deoxycytidine) is a purine nucleoside analog with potential antitumor activity for the study of apoptosis.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Solubility Information

Solubility	DMSO: 50 mg/mL (207.26 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1452 mL	20.7262 mL	41.4525 mL
5 mM	0.829 mL	4.1452 mL	8.2905 mL
10 mM	0.4145 mL	2.0726 mL	4.1452 mL
50 mM	0.0829 mL	0.4145 mL	0.829 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Robak T, Robak P. Purine nucleoside analogs in the treatment of rarer chronic lymphoid leukemias. Curr Pharm Des. 2012;18(23):3373-88.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481