


5'-DMTr-dG(iBu) -Methyl phosphonamidite

Chemical Properties

CAS No. :	115131-08-3
Formula:	C42H53N6O7P
Molecular Weight:	784.88
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside phosphonamidite
-------------	----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2741 mL	6.3704 mL	12.7408 mL
5 mM	0.2548 mL	1.2741 mL	2.5482 mL
10 mM	0.1274 mL	0.637 mL	1.2741 mL
50 mM	0.0255 mL	0.1274 mL	0.2548 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481