


DMTr-FNA-C(Bz)phosphoramidite

Chemical Properties

CAS No. : 326802-62-4
Formula: C₄₅H₅₂N₅O₈P
Molecular Weight: 821.9
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside phosphoramidites; Nucleoside Derivatives - Acyclic nucleosides;
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2167 mL	6.0835 mL	12.1669 mL
5 mM	0.2433 mL	1.2167 mL	2.4334 mL
10 mM	0.1217 mL	0.6083 mL	1.2167 mL
50 mM	0.0243 mL	0.1217 mL	0.2433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481