


N6-Ethyl-2'-deoxyadenosine

Chemical Properties

CAS No. : 136050-93-6
Formula: C₁₂H₁₇N₅O₃
Molecular Weight: 279.3
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Nucleoside Derivatives - 6-Modified purine nucleosides,
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5804 mL	17.9019 mL	35.8038 mL
5 mM	0.7161 mL	3.5804 mL	7.1608 mL
10 mM	0.358 mL	1.7902 mL	3.5804 mL
50 mM	0.0716 mL	0.358 mL	0.7161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481