


[bAla8]-Neurokinin A(4-10)

Chemical Properties

CAS No. : 122063-01-8
Formula: C₃₅H₅₆N₈O₁₀
Molecular Weight: 780.94
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	[bAla8]-Neurokinin A(4-10) is a neurokinin 2 (NK2) receptor agonist.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2805 mL	6.4025 mL	12.8051 mL
5 mM	0.2561 mL	1.2805 mL	2.561 mL
10 mM	0.1281 mL	0.6403 mL	1.2805 mL
50 mM	0.0256 mL	0.1281 mL	0.2561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481