

Secretin (28-54), human TFA

Chemical Properties

CAS No. :

Formula: C132H221N44F3O42

Molecular Weight: 3153.48

Appearance: no data available

Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Secretin (28-54), human TFA is a 27-amino acid peptide that works on the human Secretin receptor.
-------------	---

Solubility Information

Solubility	H2O: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3171 mL	1.5855 mL	3.1711 mL
5 mM	0.0634 mL	0.3171 mL	0.6342 mL
10 mM	0.0317 mL	0.1586 mL	0.3171 mL
50 mM	0.0063 mL	0.0317 mL	0.0634 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Pang RT, et al. Role of N-linked glycosylation on the function and expression of the human secretin receptor. Endocrinology. 1999 Nov;140(11):5102-11.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481