


Data Sheet (Cat.No.TP1381L)

HAEGT TFA(852155-81-8 free base)

Chemical Properties

CAS No. : TP1381L
Formula: C₂₂H₃₂F₃N₇O₁₁
Molecular Weight: 627.53
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	HAEGT TFA is the first N-terminal 1-5 residues of GLP-1 peptide.
Targets(IC ₅₀)	Others

Solubility Information

Solubility	DMSO: 10 mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5935 mL	7.9677 mL	15.9355 mL
5 mM	0.3187 mL	1.5935 mL	3.1871 mL
10 mM	0.1594 mL	0.7968 mL	1.5935 mL
50 mM	0.0319 mL	0.1594 mL	0.3187 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481