


Angiotensin I/II (1-6)

Chemical Properties

CAS No. : 47896-63-9
Formula: C₃₆H₅₅N₁₁O₁₀
Molecular Weight: 801.89
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Angiotensin I/II (1-6) is a peptide (ASP-ARG-VAL-TYR-ILE-HIS) that contains the amino acids 1-6 and is converted from Angiotensin I/II peptide.
Targets(IC50)	RAAS,Others

Solubility Information

Solubility	H ₂ O: 5 mM, DMSO: Insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2471 mL	6.2353 mL	12.4705 mL
5 mM	0.2494 mL	1.2471 mL	2.4941 mL
10 mM	0.1247 mL	0.6235 mL	1.2471 mL
50 mM	0.0249 mL	0.1247 mL	0.2494 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Erdös EG, et al. Conversion of angiotensin I to angiotensin II. Am J Med. 1976 May 31;60(6):749-59.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481