


Asp-Asp-Asp-Asp-Asp acetate

Chemical Properties

CAS No. :	TP1791L
Formula:	C ₂₂ H ₃₁ N ₅ O ₁₈
Molecular Weight:	653.51
Appearance:	no data available
Storage:	keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Asp-Asp-Asp-Asp-Asp acetate is a peptide consists of 5 Asp.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Insoluble, H ₂ O: 10 mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5302 mL	7.651 mL	15.302 mL
5 mM	0.306 mL	1.5302 mL	3.0604 mL
10 mM	0.153 mL	0.7651 mL	1.5302 mL
50 mM	0.0306 mL	0.153 mL	0.306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481