

M871

Chemical Properties

CAS No. : 908844-75-7
Formula: C108H163N27O28
Molecular Weight: 2287.64
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Selective galanin GAL2 receptor antagonist (K _i values are 13.1 and 420 nM for GAL2 and GAL1 receptors respectively). Blocks the pro-nociceptive effect of GAL2 receptor agonists.
-------------	---

Solubility Information

Solubility	H ₂ O: 1 mg/mL, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4371 mL	2.1857 mL	4.3713 mL
5 mM	0.0874 mL	0.4371 mL	0.8743 mL
10 mM	0.0437 mL	0.2186 mL	0.4371 mL
50 mM	0.0087 mL	0.0437 mL	0.0874 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Jimenez-Andrade et al (2006) Activation of peripheral galanin receptors: differential effects on nociception. Pharmacol.Biochem.Behav. 85 273 PMID:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481