## Data Sheet (Cat.No.TP2189L)



LALAN -

## V5 Epitope Tag Peptide acetate

Chemi	ical Pro	operties

CAS	No	•
CAD	110.	•

Formula:	C66H112N16O22	
Molecular Weight:	1481.65	
Appearance:	no data available	
Storage:	keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year	

<b>Biological Description</b>		
Description	V5 Epitope Tag Peptide acetate is a peptide	
Targets(IC50)	Others	

Solubility Information	
Solubility	DMSO: 10 mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6749 mL	3.3746 mL	6.7492 mL
5 mM	0.135 mL	0.6749 mL	1.3498 mL
10 mM	0.0675 mL	0.3375 mL	0.6749 mL
50 mM	0.0135 mL	0.0675 mL	0.135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

## Reference

Southern J A , Young D F , Heaney F , et al. Identification of an epitope on the P and V proteins of simian virus 5 that distinguishes between two isolates with different biological characteristics[J]. Journal of General Virology,

## Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481