


Lys-γ3-MSH(human)

Chemical Properties

CAS No. : 156159-18-1
Formula: C128H193N45O39S
Molecular Weight: 3018.25
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pro-opiomelanocortin (POMC) derived peptide that stimulates lipolysis
Targets(IC50)	Others

Solubility Information

Solubility	10% acetonitrile / H2O: 1 mg/ml, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3313 mL	1.6566 mL	3.3132 mL
5 mM	0.0663 mL	0.3313 mL	0.6626 mL
10 mM	0.0331 mL	0.1657 mL	0.3313 mL
50 mM	0.0066 mL	0.0331 mL	0.0663 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Wu W, Bi C, Bence A K, et al. Antitumor activity of Chk1 inhibitor LY2606368 as a single agent in SW1990 human pancreas orthotopic tumor model. Cancer Research, 2012, 72(8 Supplement): 1776.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481