


Leucopyrokinin

Chemical Properties

CAS No. : 104052-00-8
Formula: C42H66N12O12
Molecular Weight: 931.05
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Leucopyrokinin is a myotropic neuropeptide.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0741 mL	5.3703 mL	10.7406 mL
5 mM	0.2148 mL	1.0741 mL	2.1481 mL
10 mM	0.1074 mL	0.537 mL	1.0741 mL
50 mM	0.0215 mL	0.1074 mL	0.2148 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ryszka F, Dolińska B, Suszka-Świtek A, Rykaczewska-Czerwińska M, Konopińska D, Kuczer M, Plech A. Distribution in Rats Internal Organs of Intraperitoneally Given 125I-Labeled Heptapeptide [2-8]-Leucopyrokinin ([2-8]-LPK), a

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481