

WT-1 A1

Chemical Properties


CAS No. : 852243-39-1

Formula: C55H74N10O13S

Molecular Weight: 1115.3

Appearance: no data available

Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	WT-1 A1 is an attractive target for immunotherapy in patients with pancreatic adenocarcinoma.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8966 mL	4.4831 mL	8.9662 mL
5 mM	0.1793 mL	0.8966 mL	1.7932 mL
10 mM	0.0897 mL	0.4483 mL	0.8966 mL
50 mM	0.0179 mL	0.0897 mL	0.1793 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Terashima T, Mizukoshi E, Arai K, Yamashita T, Yoshida M, Ota H, Onishi I, Kayahara M, Ohtsubo K, Kagaya T, Honda M, Kaneko S. P53, hTERT, WT-1, and VEGFR2 are the most suitable targets for cancer vaccine therapy in

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481