


Clerodin

Chemical Properties

CAS No. : TSP-21087008
Formula: C₂₄H₃₄O₇
Molecular Weight: 434.53
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Clerodin is a useful organic compound for research related to life sciences and the catalog number is TSP-21087008.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3013 mL	11.5067 mL	23.0134 mL
5 mM	0.4603 mL	2.3013 mL	4.6027 mL
10 mM	0.2301 mL	1.1507 mL	2.3013 mL
50 mM	0.046 mL	0.2301 mL	0.4603 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481