


Injection Date : 24-May-18, 09:36:27
Sample Name : 130014
Acq. Operator : 1
Spec. Reported : Integration
Method Info :

Tgt Mass(EZX) :
Location : Vial 3
Inj : 0
Inj Volume : 5 ul

Column: Agilent C18, 100mm*4.6mm, 3.0um
Temperature: 40. Flow rate: 1.5mL/min
Mobile phase: Water: ACN 70: 30



Integration Results for DAD1 B, Sig=300,4 Ref=360,100

RetTim	Width	Area	Height	Area%
10.36	0.08	6917.66	2138.45	98.34
11.02	0.05	117.03	24.68	1.66