

Datasheet

Inhibitors / Agonists / Screening Libraries

A DRUG SCREENING EXPERT

Product Name : LY-411575 isomer 2
Catalog Number : T11904
CAS Number :
Molecular Formula : C₂₆H₂₃F₂N₃O₄
Molecular Weight : 479.48

Description: LY-411575 isomer 2 is a potent γ -secretase inhibitor. LY-411575 isomer 2 is an isomer of LY411575.

Storage: 2 years -80°C in solvent; 3 years -20°C powder;

Receptor (IC50)

Others

Reference

1. Fauq AH, et al. A multigram chemical synthesis of the gamma-secretase inhibitor LY411575 and its diastereoisomers. Bioorg Med Chem Lett. 2007 Nov 15;17(22):6392-5. Epub 2007 Aug 30.

FOR RESEARCH PURPOSES ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

Information for product storage and handling is indicated on the product datasheet. Targetmol products are stable for long term under the recommended storage conditions. Our products may be shipped under different conditions as many of them are stable in the short-term at higher or even room temperatures. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, please follow the storage recommendations on the product data sheet.