

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 6.82 - 6.89 (m, 1 H) 6.91 - 6.99 (m, 1 H) 7.06 (ddd, $J=7.84, 6.08, 2.30$ Hz, 1 H) 7.11 - 8.80 (m, 4 H) 7.38 - 7.48 (m, 2 H) 8.24 (d, $J=7.67$ Hz, 1 H) 8.48 (dd, $J=11.18, 2.63$ Hz, 1 H) 10.74 (s, 1 H) 11.81 (s, 1 H) 13.45 - 13.83 (m, 1 H)

