

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 1.45 (d, $J=6.94$ Hz, 3 H) 3.88 (s, 3 H) 5.03 (t, $J=7.26$ Hz, 1 H) 6.20 (s, 1 H) 6.69 (s, 1 H) 6.98 (d, $J=8.09$ Hz, 1 H) 7.11 (t, $J=8.86$ Hz, 2 H) 7.45 (dd, $J=8.61, 5.65$ Hz, 2 H) 7.69 (s, 1 H) 8.00 (s, 1 H) 8.03 (d, $J=2.57$ Hz, 1 H) 8.17 (dd, $J=2.57, 1.54$ Hz, 1 H) 9.22 (s, 1 H) 9.42 (s, 1 H)

